

Product Datasheet

TET2 Antibody (orb1238177)

Description	TET2 Antibody	250-
Species/Host	Rabbit	130-
Reactivity	Human, Mouse	100-
Conjugation	Unconjugated	70-
Tested Applications	ELISA, WB	Western blot analysis of TET2 in 3T3 cel...
Immunogen	TET2 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human TET2. The immunogen is located within amino acids 1550 - 1600 of TET2.	
Target	TET2	
Preservatives	TET2 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	TET2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.	
Note	For research use only	
Application notes	TET2 antibody can be used for detection of TET2 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 220 kDa Observed: 220 kDa	
Uniprot ID	Q6N021	
NCBI	NP_001120680	
Dilution Range	TET2 antibody can be used for detection of TET2 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	