

Product Datasheet

TCIRG1 Antibody (orb1238180)

Description

TCIRG1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

TCIRG1 antibody was raised against a 19 amino acid synthetic peptide near the amino terminus of human TCIRG1. The immunogen is located within amino acids 70 - 120 of TCIRG1.

Target

TCIRG1

Preservatives

TCIRG1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TCIRG1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

TCIRG1 antibody can be used for detection of TCIRG1 by Western blot at 0.5 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

91 kDa

Uniprot ID

[Q13488](#)

NCBI


[NP_006010](#)

Dilution Range

TCIRG1 antibody can be used for detection of TCIRG1 by Western blot at 0.5 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of TCIRG1 in EL4 c...