

## Product Datasheet

### SLAIN1 Antibody (orb1238210)

## Description

## SLAIN1 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

ELISA, IHC-P

### Applications

### Immunogen

SLAIN1 antibody was raised against a 16 amino acid synthetic peptide near the center of human SLAIN1. The immunogen is located within amino acids 370 - 420 of SLAIN1.

### Target

SLAIN1

### Preservatives

SLAIN1 Antibody is supplied in PBS containing 0.02% sodium azide.

### Form/Appearance

Liquid

### Concentration

1 mg/mL

### Storage

SLAIN1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

SLAIN1 antibody can be used for detection of SLAIN1 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.

### Isotype

IgG

### Clonality

Polyclonal

### Uniprot ID

[Q8ND83](#)

### NCBI


[NP\\_002456](#)

### Dilution Range

SLAIN1 antibody can be used for detection of SLAIN1 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.

### Expiration Date

12 months from date of receipt.



Immunohistochemistry of SLAIN1 in human ...