


Product Datasheet

SKA3 Antibody (orb1238211)

Description	SKA3 Antibody	A B
Species/Host	Rabbit	 Western blot analysis of SKA3 in human t...
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	SKA3 antibody was raised against a 13 amino acid synthetic peptide from near the carboxy terminus of human SKA3. The immunogen is located within amino acids 290 - 340 of SKA3.	
Target	SKA3	
Preservatives	SKA3 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	SKA3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	SKA3 antibody can be used for detection of SKA3 by Western blot at 0.5 - 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
Uniprot ID	Q8IX90	
NCBI	NP_659498	
Dilution Range	SKA3 antibody can be used for detection of SKA3 by Western blot at 0.5 - 1 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	