


Product Datasheet

SARS-CoV Spike Antibody (orb1238233)

Description	SARS-CoV Spike Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Anti-SARS-CoV Spike antibody (orb1238233) was raised against a peptide corresponding to 17 amino acids near the center of SARS-CoV Spike glycoprotein. The immunogen is located within amino acids 550-600 of SARS-CoV Spike.
Target	S
Preservatives	SARS-CoV Spike Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SARS-CoV Spike antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P59594
NCBI	P59594
Dilution Range	SARS-CoV Spike antibody can be used for the detection of SARS-CoV Spike protein in ELISA. It will detect 2 ng of free peptide at 1 µg/mL.
Expiration Date	12 months from date of receipt.



ELISA Test.
Antibodies:
SARS-CoV
Spike A...