

Product Datasheet

SARS-CoV Spike Antibody (orb1238235)

Description

SARS-CoV Spike Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

Anti-SARS-CoV Spike antibody (orb1238235) was raised against a peptide corresponding to 14 amino acids near the amino of SARS-CoV Spike glycoprotein. The immunogen is located within the first 50 amino acids of SARS-CoV Spike.

Target

S

Preservatives

SARS-CoV Spike Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV Spike antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[P59594](#)

NCBI


[P59594](#)

Dilution Range

SARS-CoV Spike antibody can be used for the detection of SARS-CoV Spike protein in ELISA. It will detect 5 ng of free peptide at 1 µg/mL.

Expiration Date

12 months from date of receipt.



ELISA Test.
Antibodies:
SARS-CoV
Spike A...