


Product Datasheet

ORF10 Antibody (orb1238240)

Description	ORF10 Antibody	 ELISA Validation. Antibodies: SARS-CoV-2...
Species/Host	Rabbit	
Reactivity	Virus	
Conjugation	Unconjugated	
Tested Applications	ELISA	
Immunogen	Anti-SARS-CoV-2 (COVID-19) ORF10 antibody (orb1238240) was raised against a peptide corresponding to 13 amino acids near the carboxyl terminus of SARS-CoV-2 (COVID-19) ORF10 protein. The immunogen is located within the last 50 amino acids of the SARS-CoV-2 (COVID-19) ORF10 protein.	
Target	ORF10	
Preservatives	SARS-CoV-2 (COVID-19) ORF10 antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	SARS-CoV-2 (COVID-19) ORF10 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	4kD	
NCBI	QSJ03260.1	
Dilution Range	Antibody validated: SARS-CoV-2 (COVID-19) ORF10 antibody can detect 2 ng of free peptide at 1 µg/mL in ELISA. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	