

Product Datasheet

Non-structural protein 6 Antibody (orb1238241)

Description

Non-structural protein 6 Antibody

Species/Host

Rabbit

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA

Applications

Immunogen

Anti-SARS-CoV-2 (COVID-19) NSP6 antibody (orb1238241) was raised against a peptide corresponding to 18 amino acids in the center of SARS-CoV-2 (COVID-19) NSP6 protein. The immunogen is located in between 230-280 amino acids of the SARS-CoV-2 (COVID-19) NSP6 protein.

Target

Non-structural protein 6

Preservatives

SARS-CoV-2 (COVID-19) NSP6 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SARS-CoV-2 (COVID-19) NSP6 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[P0DTC1](#)

NCBI


[YP_009742613.1](#)

Dilution Range

Antibody validated: SARS-CoV-2 (COVID-19) NSP6 antibody can detect 2ng of free peptide at 1 µg/mL in ELISA.

Expiration Date

12 months from date of receipt.



ELISA
Validation.
Antibodies:
SARS-CoV-
2...