

## Product Datasheet

### PTCHD2 Antibody (orb1238274)

**Description**

PTCHD2 Antibody


**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

ELISA, IHC-P

**Applications**

Immunohistochemistry  
of PTCHD2 in human

**Immunogen**

PTCHD2 antibody was raised against a 19 amino acid synthetic peptide from near the amino terminus of human PTCHD2. The immunogen is located within amino acids 210 - 260 of PTCHD2.

...

**Target**

PTCHD2

**Preservatives**

PTCHD2 Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

PTCHD2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

PTCHD2 antibody can be used for detection of PTCHD2 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.

**Isotype**

IgG

**Clonality**

Polyclonal

**Uniprot ID**
[Q9P2K9](#)
**NCBI**
[Q9P2K9](#)
**Dilution Range**

PTCHD2 antibody can be used for detection of PTCHD2 by immunohistochemistry at 5 µg/mL. Antibody validated: Immunohistochemistry in human samples. All other applications and species not yet tested.

**Expiration Date**

12 months from date of receipt.