

Product Datasheet

PST1 Antibody (orb1238276)

Description

PST1 Antibody

250-

Species/Host

Rabbit

95-

Pst1

Reactivity

Yeast

70-

Conjugation

Unconjugated

35-

Tested

ELISA, WB

Applications

Western blot
analysis of
TM yeast
Pst1 p...

Immunogen

Pst1 antibody was raised against a synthetic peptide from the amino terminus of yeast Pst1. The immunogen is located within amino acids 80 - 130 of Pst1.

Target

PST1

Preservatives

Pst1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Pst1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Pst1 antibody can be used for detection of Pst1 by Western blot and ELISA at 1 µg/mL.

Isotype

IgG

Clonality

Polyclonal

MW

49 kDa

Uniprot ID

[Q12355](#)

NCBI

[Q12355](#)

Dilution Range

Pst1 antibody can be used for detection of Pst1 by Western blot and ELISA at 1 µg/mL.

Expiration Date

12 months from date of receipt.