

Product Datasheet

PRRX1 Antibody (orb1238278)

Description

PRRX1 Antibody

31-
PRRX1

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Western blot
analysis of
PRRX1 in
HeLa c...

Immunogen

PRRX1 antibody was raised against a 19 amino acid peptide near the center of human PRRX1. The immunogen is located within amino acids 70 - 120 of PRRX1.

Target

PRRX1

Preservatives

PRRX1 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PRRX1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

PRRX1 antibody can be used for detection of PRRX1 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 27 kDa Observed: 26 kDa

Uniprot ID

[P54821](#)

NCBI

[NP_073207](#)

Dilution Range

PRRX1 antibody can be used for detection of PRRX1 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.