

Product Datasheet

PPFIA3 Antibody (orb1238279)

Description

PPFIA3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

PPFIA3 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human PPFIA3. The immunogen is located within amino acids 1100 - 1150 of PPFIA3.

Target

PPFIA3

Preservatives

PPFIA3 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PPFIA3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

PPFIA3 antibody can be used for detection of PPFIA3 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 131 kDa Observed: 131 kDa

Uniprot ID

[O75145](#)

NCBI

[NP_003651](#)

Dilution Range

PPFIA3 antibody can be used for detection of PPFIA3 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.

