

Product Datasheet

ORAI3 Antibody (orb1238293)

Description	ORAI3 Antibody	35-
Species/Host	Mouse	28-
Reactivity	Human, Rat	
Conjugation	Unconjugated	21-
Tested Applications	ELISA, WB	Western blot analysis of ORAI3 in rat lu...
Immunogen	A 19 amino acid synthetic peptide from near the carboxy terminus of human ORAI3.	
Target	ORAI3	
Preservatives	ORAI3 Monoclonal Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	ORAI3 monoclonal antibody can be stored at -20°C, stable for one year.	
Note	For research use only	
Application notes	ORAI3 antibody can be used for detection of ORAI3 by Western blot at 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.	
Isotype	IgG2b	
Clonality	Monoclonal	
Uniprot ID	Q9BRQ5	
NCBI	Q9BRQ5	
Dilution Range	ORAI3 antibody can be used for detection of ORAI3 by Western blot at 2 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	