

Product Datasheet

NRAS Antibody (orb1238297)

Description

NRAS Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

N-RAS antibody was raised against a 15 amino acid peptide near the carboxy terminus of human N-RAS. The immunogen is located within the last 50 amino acids of N-RAS.

Target

NRAS

Preservatives

N-RAS Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

N-RAS antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

N-RAS antibody can be used for detection of N-RAS by Western blot at 1 and 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 22 kDa Observed: 22 kDa

Uniprot ID

[P01111](#)

NCBI

[NP_002515](#)

Dilution Range

N-RAS antibody can be used for detection of N-RAS by Western blot at 1 and 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.

