

Product Datasheet

NOL3 Antibody (orb1238298)

Description

NOL3 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Nop30 antibody was raised against a 14 amino acid synthetic peptide from near the carboxy terminus of human Nop30. The immunogen is located within the last 50 amino acids of Nop30.

Target

NOL3

Preservatives

Nop30 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Nop30 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Nop30 antibody can be used for detection of Nop30 by Western blot at 0.5 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 0.5 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[O60936](#)

NCBI


[NP_001171986](#)

Dilution Range

Nop30 antibody can be used for detection of Nop30 by Western blot at 0.5 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 0.5 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of Nop30 in mouse ...