


Product Datasheet

MFSD6 antibody (orb12383)

Description	Goat polyclonal to MFSD6.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MFSD6
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 90kDa band observed in Human Heart lysates (calculated MW of 88.1kDa according to NP_060164.3). Recommended concentration: 1-3µg/ml.
Protein Sequence	RPTTHPTNASHQ
Clonality	Polyclonal
MW	88.1
Entrez	54842
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...