

Product Datasheet

NLRP14 Antibody (orb1238306)

Description

NLRP14 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

NALP14 antibody was raised against a 16 amino acid synthetic peptide near the amino terminus of human NALP14. The immunogen is located within amino acids 200 - 250 of NALP14.

Target

NLRP14

Preservatives

NALP14 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

NALP14 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

NALP14 antibody can be used for detection of NALP14 by Western blot at 1 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[Q86W24](#)

NCBI

[NP_789792](#)

Dilution Range

NALP14 antibody can be used for detection of NALP14 by Western blot at 1 µg/mL. Antibody validated: Western Blot in rat samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.

