


Product Datasheet

MUC1 Antibody (orb1238310)

Description	MUC1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	MUC1 antibody was raised against a 17 amino acid peptide near the carboxy terminus of human MUC1. The immunogen is located within the last 50 amino acids of MUC1 .	Western blot analysis of MUC1 in HeLa ce...
Target	MUC1	
Preservatives	MUC1 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	MUC1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.	
Note	For research use only	
Application notes	MUC1 antibody can be used for detection of MUC1 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 53 (+ multiple others) kDa Observed: 42, 120 kDa	
Uniprot ID	P15941	
NCBI	NP_001191215	
Dilution Range	MUC1 antibody can be used for detection of MUC1 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	