

Product Datasheet

NKIRAS2 Antibody (orb1238342)

Description

NKIRAS2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

KappaB ras2 antibody was raised against a 14 amino acid synthetic peptide from near the carboxy terminus of human KappaB ras2. The immunogen is located within the last 50 amino acids of KappaB ras2.

Target

NKIRAS2

Preservatives

KappaB ras2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

KappaB ras2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

KappaB ras2 antibody can be used for detection of KappaB ras1 by Western blot at 1 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID

[Q9NYR9](#)

NCBI


[NP_060065](#)

Dilution Range

KappaB ras2 antibody can be used for detection of KappaB ras1 by Western blot at 1 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of KappaB ras2 in ...