

Product Datasheet

JAKMIP2 Antibody (orb1238345)

Description

JAKMIP2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Rabbit polyclonal JAKMIP2 antibody was raised against a 16 amino acid peptide near the center of human JAKMIP2. The immunogen is located within amino acids 250 - 300 of JAKMIP2.

Target

JAKMIP2

Preservatives

JAKMIP2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

JAKMIP2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

JAKMIP2 antibody can be used for detection of JAKMIP2 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 89 kDa

Uniprot ID

[Q96AA8](#)

NCBI


[NP_055605](#)

Dilution Range

JAKMIP2 antibody can be used for detection of JAKMIP2 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of JAKMIP2 in mouse...