


Product Datasheet

ITCH Antibody (orb1238347)

Description	ITCH Antibody	
Species/Host	Rabbit	 <p>Western blot analysis of ITCH in 3T3 cel...</p>
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	ITCH antibody was raised against a 15 amino acid peptide near the amino terminus of human ITCH .The immunogen is located within amino acids 40 - 90 of ITCH.	
Target	ITCH	
Preservatives	ITCH antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	ITCH antibody can be stored at 4°C for three months and -20°C, stable for up to one year.	
Note	For research use only	
Application notes	ITCH antibody can be used for detection of ITCH by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 99 kDa Observed: 99 kDa	
Uniprot ID	Q96J02	
NCBI	NP_001244066	
Dilution Range	ITCH antibody can be used for detection of ITCH by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	