

Product Datasheet

IL17RE Antibody (orb1238355)

Description

IL17RE Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Rabbit polyclonal IL-17RE antibody was raised against a 16 amino acid peptide near the center of human IL-17RE. The immunogen is located within amino acids 230 - 280 of IL-17RE.

Target

IL17RE

Preservatives

IL-17RE Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

IL-17RE antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

IL-17RE antibody can be used for detection of IL-17RE by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 78 kDa

Uniprot ID

[Q8NFR9](#)

NCBI


[NP_705616](#)

Dilution Range

IL-17RE antibody can be used for detection of IL-17RE by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of IL-17RE in human samples.