

Product Datasheet

HEPACAM2 Antibody (orb1238367)

Description

HEPACAM2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Rabbit polyclonal HEPACAM2 antibody was raised against a 19 amino acid peptide near the carboxy terminus of human HEPACAM2. The immunogen is located within amino acids 370 - 420 of HEPACAM2.

Target

HEPACAM2

Preservatives

HEPACAM2 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

HEPACAM2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

HEPACAM2 antibody can be used for detection of HEPACAM2 by Western blot at 0.5 - 1 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

51 kDa

Uniprot ID

[A8MVW5](#)

NCBI


[NP_001034461](#)

Dilution Range

HEPACAM2 antibody can be used for detection of HEPACAM2 by Western blot at 0.5 - 1 µg/mL. Antibody validated: Western Blot in mouse samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of HEPACAM2 in mouse samples.