

Product Datasheet

TRIM6 antibody (orb12384)

Description	Goat polyclonal to TRIM6.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	TRIM6 (aa118-131)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize reported isoforms 1 and 2 NP_001003818.1;NP_477514.1. Reported variants represent identical protein: NP_001185573.1, NP_001185574.1
Protein Sequence	KQLKAVLCADHGEK
Clonality	Polyclonal
MW	59.4; 56.4; 36.4
Entrez	117854 , 445372
Dilution Range	ELISA: 1:8000, WB: 0.2 µg/ml
Expiration Date	12 months from date of receipt.