


## Product Datasheet

### DNAJA4 Antibody (orb1238422)

| Description         | DNAJA4 Antibody                                                                                                                                                                           |  |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Species/Host        | Rabbit                                                                                                                                                                                    |                                                                                     |
| Reactivity          | Human, Mouse, Rat                                                                                                                                                                         |                                                                                     |
| Conjugation         | Unconjugated                                                                                                                                                                              |                                                                                     |
| Tested Applications | ELISA, WB                                                                                                                                                                                 |                                                                                     |
| Immunogen           | DNAJA4 antibody was raised against a 15 amino acid peptide near the center of human DNAJA4. The immunogen is located within amino acids 210 - 260 of DNAJA4.                              | Western blot analysis of DNAJA4 in human...                                         |
| Target              | DNAJA4                                                                                                                                                                                    |                                                                                     |
| Preservatives       | DNAJA4 antibody is supplied in PBS containing 0.02% sodium azide.                                                                                                                         |                                                                                     |
| Form/Appearance     | Liquid                                                                                                                                                                                    |                                                                                     |
| Concentration       | 1 mg/mL                                                                                                                                                                                   |                                                                                     |
| Storage             | DNAJA4 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.                                                                                               |                                                                                     |
| Note                | For research use only                                                                                                                                                                     |                                                                                     |
| Application notes   | DNAJA4 antibody can be used for detection of DNAJA4 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested. |                                                                                     |
| Isotype             | IgG                                                                                                                                                                                       |                                                                                     |
| Clonality           | Polyclonal                                                                                                                                                                                |                                                                                     |
| MW                  | Predicted: 41, 44, 47 kDa Observed: 40 kDa                                                                                                                                                |                                                                                     |
| Uniprot ID          | <a href="#">Q8WW22</a>                                                                                                                                                                    |                                                                                     |
| NCBI                | <a href="#">NP_061072</a>                                                                                                                                                                 |                                                                                     |
| Dilution Range      | DNAJA4 antibody can be used for detection of DNAJA4 by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human samples. All other applications and species not yet tested. |                                                                                     |
| Expiration Date     | 12 months from date of receipt.                                                                                                                                                           |                                                                                     |