


## Product Datasheet

### BID Antibody (orb1238472)

Description	BID Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	ELISA, WB	
Immunogen	Bid antibody was raised against a peptide corresponding to 14 amino acids near the middle of human Bid. The immunogen is located within amino acids 100 - 150 of Bid.	Western blot analysis of Bid in A549 cel...
Target	BID	
Preservatives	Bid Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	Bid antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	Bid antibody can be used for detection of Bid by Western blot at 0.5 to 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Isotype	IgG	
Clonality	Polyclonal	
Uniprot ID	<a href="#">P55957</a>	
NCBI	<a href="#">AAH36364</a>	
Dilution Range	Bid antibody can be used for detection of Bid by Western blot at 0.5 to 2 µg/mL. Antibody validated: Western Blot in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	