

Product Datasheet

TUBA1A Antibody (Biotin) (orb1238506)

Description

TUBA1A Antibody (Biotin)

Species/Host

Rabbit

Reactivity

Gallus, Human, Mouse, Rabbit, Rat, Zebrafish

Conjugation

Biotin

Tested

ELISA, WB

Applications

Immunogen

Biotin-Alpha-tubulin antibody was raised against an 18 amino acid peptide near the carboxy terminus of human alpha-tubulin

Target

TUBA1A

Preservatives

Biotin-Alpha-tubulin antibody is supplied in PBS containing 1% rAlbumin and 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Biotin-Alpha-tubulin antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

Biotin-Alpha-tubulin antibody can be used for detection of alpha-tubulin by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human, mouse, and rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 50 kDa Observed: 55 kDa

Uniprot ID

[Q71U36](#)

NCBI


[CAA25855](#)

Dilution Range

Biotin-Alpha-tubulin antibody can be used for detection of alpha-tubulin by Western blot at 1 - 2 µg/ml. Antibody validated: Western Blot in human, mouse, and rat samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of Alpha-tubulin i...