

Product Datasheet

ZNF821 Antibody (orb1238524)

Description

ZNF821 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, WB

Applications
Immunogen

ZNF821 antibody was raised against an 18 amino acid peptide near the carboxy terminus of human ZNF821. The immunogen is located within amino acids 270 - 320 of ZNF821.

Target

ZNF821

Preservatives

ZNF821 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

ZNF821 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

ZNF821 antibody can be used for detection of ZNF821 by Western blot at 0.5 - 1 µg/ml. Antibody can also be used for Immunocytochemistry at 2.5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

MW

Predicted: 41 kDa Observed: 38 kDa

Uniprot ID
[O75541](#)
NCBI
[NP_001188482.1](#)
Dilution Range


ZNF821 antibody can be used for detection of ZNF821 by Western blot at 0.5 - 1 µg/ml. Antibody can also be used for Immunocytochemistry at 2.5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.



Western blot analysis of ZNF821 in Jurka...



Immunocytochemistry of ZNF821 in Jurkat ...



Immunofluorescence of ZNF821 in Jurkat ce...