

Product Datasheet

ZBTB7A Antibody (orb1238561)

Description

ZBTB7A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

ZBTB7A antibody was raised against a 16 amino acid synthetic peptide near the amino terminus of human ZBTB7A. The immunogen is located within the first 50 amino acids of ZBTB7A.

Target

ZBTB7A

Preservatives

ZBTB7A Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

ZBTB7A antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

ZBTB7A antibody can be used for detection of ZBTB7A by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2.5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

Uniprot ID


[O95365](#)

NCBI


[AAH84568](#)

Dilution Range

ZBTB7A antibody can be used for detection of ZBTB7A by Western blot at 1 - 2 µg/mL. Antibody can also be used for immunohistochemistry starting at 2.5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples




Western blot analysis of ZBTB7A in human...



Immunohistochemistry of ZBTB7A in human

...



Immunofluorescence of ZBTB7A in human ov...