

Product Datasheet

UNG Antibody (orb1238581)

Description

UNG Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested Applications

ELISA, IF, IHC-P, WB

Immunogen

Anti-UNG1 antibody (orb1238581) was raised against a peptide corresponding to 14 amino acids near the center of human UNG1. The immunogen is located within amino acids 200-250 of UNG1.

Target

UNG

Preservatives

UNG1 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

UNG1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 34kD
Observed: 34 kD

Uniprot ID

[P13051](#)

NCBI

[NP_003353](#)


Dilution Range

WB: 0.5-4 µg/mL; IHC: 2 µg/mL; IF: 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.


Expiration Date

12 months from date of receipt.


A B C



Western Blot
Validation in Mouse
C2C12 C...



Independent Antibody
Validation (IAV) vi...



Immunohistochemistry
Validation of UNG1 ...