

Product Datasheet

WNT5B Antibody (orb1238587)

Description

WNT5B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

WNT5 antibody was raised against a 15 amino acid peptide near the center of human WNT5B. The immunogen is located within amino acids 130 - 180 of WNT5.

Target

WNT5B

Preservatives

WNT5 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

WNT5 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

WNT5 antibody can be used for detection of WNT5 by Western blot at 1 - 2 µg/ml. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in rat samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 39 kDa Observed: 465 kDa

Uniprot ID

[Q9H1J7](#)

NCBI


[NP_110402](#)

Dilution Range


WNT5 antibody can be used for detection of WNT5 by Western blot at 1 - 2 µg/ml. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples and Immunohistochemistry in rat samples. All other applications and species not yet tested.

Expiration Date

12 months from date of receipt.



Western blot analysis of WNT5 in Jurkat ...



Immunohistochemistry of WNT5 in rat colo...