


Product Datasheet

VTI1b Antibody (orb1238600)

Description		VTI1b Antibody		<div><div>31-</div><div>A</div><div>B</div></div>	
Species/Host	Rabbit	<div><div>24-</div><div></div><div>14-</div><div></div><div>8-</div><div></div></div>			
Reactivity	Human, Mouse, Rat				
Conjugation	Unconjugated				
Tested Applications	ELISA, IHC-P, WB	Western blot analysis of VTI1b in HeLa c...			
Immunogen	VTI1b antibody was raised against a 17 amino acid peptide near the center of human VTI1b.The immunogen is located within amino acids 150 - 200 of VTI1b.				
Target	VTI1b				
Preservatives	VTI1b antibody is supplied in PBS containing 0.02% sodium azide.				
Form/Appearance	Liquid	Immunohistochemistry of VTI1b in human b...			
Concentration	1 mg/mL				
Storage	VTI1b antibody can be stored at 4°C for three months and -20°C, stable for up to one year.				
Note	For research use only				
Application notes	VTI1b antibody can be used for detection of VTI1b by Western blot at 1 - 2 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL.Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested.				
Isotype	IgG				
Clonality	Polyclonal				
MW	Predicted: 14, 19, 26 kDa Observed: 27 kDa				
Uniprot ID	Q9UEU0				
NCBI	NP_006361				
Dilution Range	VTI1b antibody can be used for detection of VTI1b by Western blot at 1 - 2 µg/ml. Antibody can also be used for immunohistochemistry starting at 5 µg/mL.Antibody validated: Western Blot in human samples and Immunohistochemistry in human samples. All other applications and species not yet tested.				
Expiration Date	12 months from date of receipt.				