



Product Datasheet

TRIM21 Antibody (orb1238609)

Description	TRIM21 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IF, IHC-P, WB	
Immunogen	Anti-TRIM21 antibody (orb1238609) was raised against a peptide corresponding to 16 amino acids near the center of human TRIM21. The immunogen is located within amino acids 290-340 of TRIM21.	 Western Blot Validation in Human Lung Ti...
Target	TRIM21	
Preservatives	TRIM21 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	TRIM21 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.	
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 54kDObserved: 54 kD	
Uniprot ID	P19474	
NCBI	NP_003132	
Dilution Range	WB: 0.5-2 µg/mL; IHC: 5 µg/mL; IF: 20 µg/mL.Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples; Immunofluorescence in human samples. All other applications and species not yet tested.	 Independent Antibody Validation (IAV) vi... Immunohistochemistry Validation of TRIM2...
Expiration Date	12 months from date of receipt.	