


Product Datasheet

SIGLEC15 Antibody (orb1238667)


Description

SIGLEC15 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Anti-SIGLEC15 antibody (orb1238667) was raised against a peptide corresponding to 14 amino acids near the carboxy terminus of human SIGLEC15. The immunogen is located within amino acids 270-320 of SIGLEC15.
Target	SIGLEC15
Preservatives	SIGLEC15 Antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SIGLEC15 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Predicted: 36kD Observed: 34-36kD
Uniprot ID	Q6ZMC9
NCBI	NP_998767
Dilution Range	WB: 2 µg/mL. Antibody validated: Western Blot in human, mouse and rat samples. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



Western Blot Validation in Human and Mouse...



Western Blot Validation in Human Tissues...



Western Blot Validation in Mouse Tissues...