

Product Datasheet

TRIM33 Antibody (orb1238680)

Description

TRIM33 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

TRIM33 antibody was raised against an 18 amino acid peptide near the center of human TRIM33. The immunogen is located within amino acids 580 - 630 of TRIM33.

Target

TRIM33

Preservatives

TRIM33 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

TRIM33 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

TRIM33 antibody can be used for detection of TRIM33 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 70, 124 kDa Observed: 70 kDa

Uniprot ID

[Q9UPN9](#)

NCBI


[NP_056990](#)

Dilution Range


TRIM33 antibody can be used for detection of TRIM33 by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western blot analysis of TRIM33 in human...



Immunohistochemistry of TRIM33 in human ...



Immunofluorescence of TRIM33 in human li...