

Product Datasheet

SARS-CoV-2 (COVID-19) Spike P681R Antibody [5H4C5] (Delta Variant) (orb1238700)

Description

SARS-CoV-2 (COVID-19) Spike

	Variant)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Anti-SARS-CoV-2 Spike 681R (Delta Variant) antibody (orb1238700) was raised against a peptide corresponding to 12 amino acids near the cleavage site of SARS-CoV-2 Delta variant Spike protein. The immunogen is located within 650-700 amino acids of SARS-CoV-2 Delta variant Spike protein.
Target	S
Preservatives	SARS-CoV-2 Spike P681R (Delta Variant) antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SARS-CoV-2 Spike P681R (Delta Variant) antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P0DTC2
NCBI	QHD43416

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Delta Variant Spike S1 protein,