

## Product Datasheet

### TNFSF4 Antibody (orb1238712)

**Description**

TNFSF4 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

ELISA, IF, IHC-P, WB

**Applications**

**Immunogen**

Rabbit polyclonal TNFSF4 antibody was raised against a 14 amino acid peptide near the carboxy terminus of human TNFSF4. The immunogen is located within the last 50 amino acids of TNFSF4.

**Target**

TNFSF4

**Preservatives**

TNFSF4 Antibody is supplied in PBS containing 0.02% sodium azide.

**Form/Appearance**

Liquid

**Concentration**

1 mg/mL

**Storage**

TNFSF4 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

**Note**

For research use only

**Application notes**

TNFSF4 antibody can be used for detection of TNFSF4 by Western blot at 0.5 - 1 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

20 kDa

**Uniprot ID**

**P23510**

**NCBI**


**NP\_003317**

**Dilution Range**


TNFSF4 antibody can be used for detection of TNFSF4 by Western blot at 0.5 - 1 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

**Expiration Date**


12 months from date of receipt.



Western blot analysis of TNFSF4 in rat s...



Immunohistochemistry of TNFSF4 in human ...



Immunofluorescence of TNFSF4 in human sp...