

Product Datasheet

PPARGC1A Antibody (orb1238778)

Description

PPARGC1A Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

PPARGC1A antibody was raised against a 19 amino acid peptide near the carboxy terminus of human PPARGC1A. The immunogen is located within the last 50 amino acids of PPARGC1A.

Target

PPARGC1A

Preservatives

PPARGC1A antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PPARGC1A antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

PPARGC1A antibody can be used for detection of PPARGC1A by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 90 kDa Observed: 88 kDa

Uniprot ID


[Q9UBK2](#)

NCBI


[NP_037393](#)

Dilution Range


PPARGC1A antibody can be used for detection of PPARGC1A by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry starting at 5 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human and mouse samples; Immunohistochemistry in



Western blot analysis of PPARGC1A in mou...



Immunohistochemistry of PPARGC1A in huma...



Immunofluorescence of PPARGC1A in human ...