




Product Datasheet

Perp Antibody (orb1238785)

Description	Perp Antibody	
Species/Host	Rabbit	 Independent Antibody Validation (IAV) vi...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, ICC, WB	
Immunogen	Anti-PERP antibody (orb1238785) was raised against a peptide corresponding to 17 amino acids near the carboxy terminus of human PERP. The immunogen is located within the last 50 amino acids of PERP.	
Target	Perp	 Western Blot Validation in Human Cell Li...
Preservatives	PERP Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	PERP antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	 Western Blot Validation in Human A431 wh...
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 21kD Observed: 21 kD	
Uniprot ID	Q9JK95	
NCBI	NP_071404	
Dilution Range	WB: 1-2 µg/mL; ICC: 10 µg/mL. Antibody validated: Western Blot in human samples and Immunocytochemistry in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	