

Product Datasheet

PACS1 Antibody (orb1238788)

Description

PACS1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

PACS1 antibody was raised against a 17 amino acid peptide near the center of human PACS1. The immunogen is located within amino acids 430 - 480 of PACS1.

Target

PACS1

Preservatives

PACS1 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

PACS1 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note

For research use only

Application notes

PACS1 antibody can be used for detection of PACS1 by Western blot at 1 - 2 µg/ml. Antibody can also be used for immunohistochemistry starting at 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in mouse samples and Immunofluorescence in human and mouse samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 106 kDa Observed: 105 kDa

Uniprot ID


[Q6VY07](#)

NCBI


[NP_060496](#)

Dilution Range


PACS1 antibody can be used for detection of PACS1 by Western blot at 1 - 2 µg/ml. Antibody can also be used for immunohistochemistry starting at 2 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in rat samples; Immunohistochemistry in mouse samples and Immunofluorescence in human and mouse samples. All other applications and species not yet tested.



Western blot analysis of PACS1 in rat brain tissue.



Immunohistochemistry of PACS1 in mouse brain tissue.



Immunofluorescence of PACS1 in human brain tissue.