

Product Datasheet

TICRR Antibody (orb1238799)

Description

TICRR Antibody

250-  TICRR

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Western blot analysis
of TICRR in human ...

Immunogen

TICRR antibody was raised against a 19 amino acid peptide near the carboxy terminus of human TICRR. The immunogen is located within amino acids 1840 - 1890 of TICRR.

Target

TICRR

Preservatives

TICRR antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage


TICRR antibody can be stored at 4°C for three months and -20°C, stable for up to one year.

Note


For research use only

Application notes

TICRR antibody can be used for detection of TICRR by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.



Immunohistochemistry
of TICRR in human s...



Immunofluorescence
of TICRR in human
sma...

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 210 kDa Observed: 220 kDa

Uniprot ID

[Q7Z2Z1](#)

NCBI

[NP_689472](#)

Dilution Range

TICRR antibody can be used for detection of TICRR by Western blot at 1 - 2 µg/ml. Antibody can also be used for Immunohistochemistry at 5 µg/mL. For Immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.