




Product Datasheet

KREMEN2 Antibody (orb1238862)

Description	KREMEN2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse, Rat	
Conjugation	Unconjugated	
Tested Applications	ELISA, ICC, WB	
Immunogen	Anti-KREMEN2 antibody (orb1238862) was raised against a peptide corresponding to 14 amino acids near the carboxy terminus of human KREMEN2. The immunogen is located within amino acids 390-440 of KREMEN2.	Independent Antibody Validation (IAV) vi...
Target	KREMEN2	
Preservatives	KREMEN2 Antibody is supplied in PBS containing 0.02% sodium azide.	
Form/Appearance	Liquid	
Concentration	1 mg/mL	
Storage	KREMEN2 antibody can be stored at 4°C for three months and -20°C, stable for up to one year.	Western Blot Validation in Human Cell Li...
Note	For research use only	
Isotype	IgG	
Clonality	Polyclonal	
MW	Predicted: 49kDObserved: 49kD	
Uniprot ID	Q8NCW0	Western Blot Validation in Mouse Spleen ...
NCBI	NP_757384	
Dilution Range	WB: 1-2 µg/mL; ICC: 5 µg/mL.Antibody validated: Western Blot in human, mouse and rat samples; Immunocytochemistry in human samples. All other applications and species not yet tested.	
Expiration Date	12 months from date of receipt.	