

Product Datasheet

ABCC9 antibody (orb12389)

Description	Goat polyclonal to ABCC9.
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	ABCC9 / SUR2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize isoforms SUR2A NP_005682.2 and SUR2B NP_064693.2.
Protein Sequence	TSEYSINNTGKAD
Clonality	Polyclonal
MW	170; 174; 174
Entrez	10060
Dilution Range	ELISA: 1:32000, WB: 0.1 µg/ml
Expiration Date	12 months from date of receipt.