

Product Datasheet

IL17RA Antibody (orb1238913)

Description

IL17RA Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications

Immunogen

IL-17RA antibody was raised against an 18 amino acid peptide near the amino terminus of human IL-17RA. The immunogen is located within amino acids 150 - 200 of IL-17RA.

Target

IL17RA

Preservatives

IL-17RA Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

IL-17RA antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

IL-17RA antibody can be used for detection of IL-17RA by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 95 kDa Observed: 135 kDa

Uniprot ID

[Q96F46](#)

NCBI


[NP_055154](#)

Dilution Range


IL-17RA antibody can be used for detection of IL-17RA by Western blot at 1 - 2 µg/mL. Antibody validated: Western Blot in mouse samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western blot analysis of IL17RA in A20 c...



Immunohistochemistry of IL-17RA in human...



Immunofluorescence of IL-17RA in human S...