

## Product Datasheet

### GFRA1 Antibody (orb1238923)

## Description

## GFRA1 Antibody

### Species/Host

Rabbit

### Conjugation

Unconjugated

### Tested

### Applications

ELISA, IF, IHC-P, WB

### Immunogen

GFR alpha 1 antibody was raised against a 14 amino acid peptide near the center of human GFR alpha 1. The immunogen is located within amino acids 360 - 410 of GFR alpha 1.

### Target

GFRA1

### Preservatives

GFR alpha 1 Antibody is supplied in PBS containing 0.02% sodium azide.

### Form/Appearance

Liquid

### Concentration

1 mg/mL

### Storage

GFR alpha 1 antibody can be stored at 4°C up to one year. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

GFR alpha 1 antibody can be used for detection of GFR alpha 1 by Western blot at 1 µg/mL. Antibody can also be used for immunohistochemistry starting at 1 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

### Isotype

IgG

### Clonality

Polyclonal

### MW

51 kDa

### Uniprot ID


[P56159](#)

### NCBI


[P56159](#)

### Dilution Range


GFR alpha 1 antibody can be used for detection of GFR alpha 1 by Western blot at 1 µg/mL. Antibody can also be used for immunohistochemistry starting at 1 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.



Western blot analysis of GFR alpha 1 in ...



Immunohistochemistry of GFR alpha 1 in h...



Immunofluorescence of GFR alpha 1 in hum...