

Product Datasheet

SLCO1B1 Antibody (orb1238941)

Description

SLCO1B1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, WB

Applications

Immunogen

SLCO1B1 antibody was raised against a 14 amino acid peptide near the carboxy terminus of human SLCO1B1. The immunogen is located within the last 50 amino acids of SLCO1B1.

Target

SLCO1B1

Preservatives

SLCO1B1 antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SLCO1B1 antibody can be stored at 4°C, stable for up to one year.

Note

For research use only

Application notes

SLCO1B1 antibody can be used for detection of SLCO1B1 by Western blot at 1 - 2 µg/ml. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality

Polyclonal

MW

Predicted: 76 kDa Observed: 75 kDa

Uniprot ID

[Q9Y6L6](#)

NCBI


[NP_006437](#)

Dilution Range


SLCO1B1 antibody can be used for detection of SLCO1B1 by Western blot at 1 - 2 µg/ml. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunocytochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Expiration Date


12 months from date of receipt.



Western blot analysis of SLCO1B1 in human samples.



Immunocytochemistry of SLCO1B1 in HepG2 cells.



Immunofluorescence of SLCO1B1 in human liver tissue.