

Product Datasheet

SLC30A8 Antibody (orb1238947)

Description

SLC30A8 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IF, IHC-P, WB

Applications
Immunogen

SLC30A8 antibody was raised against a peptide corresponding to 14 amino acids near the amino terminus of human SLC30A8.

Target

SLC30A8

Preservatives

SLC30A8 Antibody is supplied in PBS containing 0.02% sodium azide.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

SLC30A8 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

SLC30A8 antibody can be used for Western blot at 1 - 2 µg/mL. Antibody can also be used for Immunohistochemistry starting at 10 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and Immunofluorescence in human samples. All other applications and species not yet tested.

Isotype

IgG

Clonality


Polyclonal

MW


Predicted: 35, 41 kDa Observed: 35, 41 kDa

Uniprot ID
[Q8IWU4](#)
NCBI
[NP_776250](#)
Dilution Range


SLC30A8 antibody can be used for Western blot at 1 - 2 µg/mL. Antibody can also be used for Immunohistochemistry starting at 10 µg/mL. For immunofluorescence start at 20 µg/mL. Antibody validated: Western Blot in human samples; Immunohistochemistry in human samples and



Western blot analysis of SLC30A8 in huma...



Immunohistochemistry of SLC30A8 in human...



Immunofluorescence of SLC30A8 in human p...