




Product Datasheet

EBI3 Antibody (orb1238976)

Description	EBI3 Antibody	KDa	1	2	3
Species/Host	Rabbit	250 -			
Reactivity	Human	95 -			
Conjugation	Unconjugated	72 -			
Tested Applications	ELISA, IF, IHC-P, WB	55 -			
Immunogen	Anti-EBI-3 antibody (orb1238976) was raised against a peptide corresponding to 16 amino acids near the carboxy terminus of mature human EBI-3. The immunogen is located within the last 50 amino acids of EBI3.	36 -			
Target	EBI3	28 -			
Preservatives	EBI3 Antibody is supplied in PBS containing 0.02% sodium azide.				
Form/Appearance	Liquid	Recombinant Protein Validation. Loading:...			
Concentration	1 mg/mL				
Storage	EBI3 antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Immunohistochemistry Validation of EBI-3...			
Note	For research use only				
Isotype	IgG	Immunofluorescence Validation of EBI-3 i...			
Clonality	Polyclonal				
MW	Predicted: 25kDObserved: 30 kD (25kD + 2 N-Glycosylations)				
Uniprot ID	Q14213				
NCBI	NP_005746				
Dilution Range	WB: 1-4 µg/mL; IF: 20 µg/mL ; IHC: 2.5 µg/mL;.Antibody validated: Western Blot, Immunofluorescence, and Immunohistochemistry in human samples. All other applications and species not yet tested.				
Expiration Date	12 months from date of receipt.				