


Product Datasheet

SARS-CoV-2 Spike P26S Antibody (Gamma Variant) (orb1239024)


Explore Bioreagents.

SARS-CoV-2 Spike P26S Antibody (Gamma


Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Anti-SARS-CoV-2 Spike P26S (Gamma Variant) antibody (orb1239024) was raised against a peptide corresponding to 11 amino acids near the amino terminus of SARS-CoV-2 Gamma variant (P.1) Spike protein. The immunogen is located within the first 50 amino acids of SARS-CoV-2 Gamma variant (P.1) Spike protein.
Target	S
Preservatives	SARS-CoV-2 Spike P26S antibody is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SARS-CoV-2 Spike P26S antibody can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0DTC2
NCBI	QHD43416
Dilution Range	WB: 1 µg/mL.Antibody validated: Western Blot in human samples. Anti-SARS-CoV-2 Spike P26S (Gamma Variant) antibody specifically detects SARS-CoV-2 Gamma Variant (P.1) Spike S1 protein, but not SARS-CoV-2 WT and other variant Spike S1 protein by ELISA. All other applications and species not yet tested.
Expiration Date	12 months from date of receipt.



SARS-Cov-
2 Spike
P26S
(Gamma
Variant)
An...



ELISA
Validation
of Gamma
Variant
Spike ...



WB
Validation
of Gamma
Variant
Spike S1 ...