


Product Datasheet

SARS-CoV-2 (COVID-19) Spike 681P Antibody (biotin) (orb1239025)


Description

SARS-CoV-2 (COVID-19) Spike 681P Antibody

Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Immunogen	Anti-SARS-CoV-2 (COVID-19) Spike Antibody (cleavage site) (biotin) (orb1239025) was raised against a peptide corresponding to 12 amino acids near the center of SARS-CoV-2 (COVID-19) Spike glycoprotein. The immunogen is located within 650-700 amino acids of SARS-CoV-2 (COVID-19) Spike protein.
Target	S
Preservatives	SARS-CoV-2 (COVID-19) Spike Antibody (cleavage site) (biotin) is supplied in PBS containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	SARS-CoV-2 (COVID-19) Spike antibody (cleavage site) (biotin) can be stored at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0DTC2
NCBI	QHD43416
Expiration Date	12 months from date of receipt.



Sandwich ELISA for SARS-CoV-2 (COVID-19)...



Sandwich ELISA for SARS-CoV-2 (COVID-19)...